




UltraStim[®]

NEUROSTIMULATION
ELECTRODES

*The Next
Generation in
Electrode Design*



PART #	SHAPE	DESCRIPTION	PACKAGES PER CASE	ELECTRODES PER CASE
USX2020		UltraStim X Electrode 2.0" x 2.0", 5cm x 5cm Square, 4 per package	10	40
USX2040		UltraStim X Electrode 2.0" x 4.0", 5cm x 10 cm Rectangle, 2 per package	10	20

To learn more, contact:

 **AXELGAARD**
MANUFACTURING CO., LTD.

HEADQUARTERS: Fallbrook, CA 92028, USA
FAX: +1 (760) 723-2356

EUROPE: DK-8520 Lystrup, Denmark
FAX: +45 8674 2281

WEBSITE: www.axelgaard.com

E-MAIL: support@axelgaard.com

UltraStim and MultiStick are Registered Trademarks
of Axelgaard Manufacturing Co., Ltd.

01/11 ©2011 Axelgaard Mfg. Co., Ltd.



ISO 13485:2003
FM 40363



New from **Axelgaard** ...

UltraStim[®]

NEUROSTIMULATION
ELECTRODES

Advanced Technology for
Better Performance and Value



*The Next Generation
in Electrode Design*

 **AXELGAARD**
MANUFACTURING CO., LTD.

The Next Generation in Electrode Design

UltraStim[®]X

NEUROSTIMULATION ELECTRODES



Design eliminates electrode stiffness for continual contact during treatment.

UltraStim X Patented Neurostimulation Electrodes –

The next generation in carbon film electrode design:

- Unique silver pattern controls current and eliminates the need for a centered lead wire.
- Made with patented **MultiStick[®] Gel**, a proven hydrogel with uncompromised adhesion and performance.
- Offset lead wire contributes to unequalled conformity in a 2" x 4" carbon film electrode.

Performance and Value – Axelgaard's advanced manufacturing techniques and superior design add flexibility, durability and economy:

- Cutting-edge processing for **UltraStim X** results in exceptional quality with exceptional value.
- The **UltraStim X** patented electrodes are available in the two essential sizes.
- Streamlined product line makes ordering simple.


The **UltraStim[®]X** Difference –

Conductive material and silver grid pattern optimize current control.

Current roll off design enhances stimulation comfort.

Silver grid pattern feature removes need for center lead wire.

Offset lead wire assures maximum conformity.



Oversize adhesive border provides edge for easy removal.

Clinician Preferred – The clinician inspired design makes **UltraStim X** the clinician's choice:

- Optimized skin adhesion during contraction for better patient comfort.
- Lead wire placement enhances skin contact during treatment by maximizing flexibility and minimizing edge lifting from stiffness.
- Oversize border to facilitate ease of proper removal, which promotes electrode durability.

UltraStim[®]X
NEUROSTIMULATION
ELECTRODES

another example of

Axelgaard's commitment to

**better results through
better technology.**

