



Ionto480[®]

Buffered
Iontophoresis Electrode Treatment Kit







Each Case Contains:

- 1 Instruction Manual
- 10 Individually Packaged Kits

Each Kit Contains:

- 1 Active Drug Delivery Electrode (size indicates fill volume)
- 1 Buffered Return Electrode
- 2 Alcohol Swabs

Part #	SHAPE	KIT CONTENTS	#KITS/ CS
DDS275		1 Small Square Active Drug Delivery Electrode (Approximate Fill Volume – 1.5cc*) 1 Buffer Electrode, 2 Alcohol Swabs	10
DDM325		1 Medium Square Active Drug Delivery Electrode (Approximate Fill Volume – 2.5cc*) 1 Buffer Electrode, 2 Alcohol Swabs	10
DDL350		1 Large Square Active Drug Delivery Electrode (Approximate Fill Volume – 4.0cc*) 1 Buffer Electrode, 2 Alcohol Swabs	10
DDB375		1 Butterfly Active Drug Delivery Electrode (Approximate Fill Volume – 2.0cc*) 1 Buffer Electrode, 2 Alcohol Swabs	10

* Amount to fully saturate the electrode may vary slightly – do not over-saturate.

To learn more, contact:



HEADQUARTERS: Fallbrook, CA 92028, USA

FAX: +1 (760) 723-2356

EUROPE: DK-8520 Lystrup, Denmark

FAX: +45 8674 2281

WEBSITE: www.axelgaard.com

E-MAIL: support@axelgaard.com



ISO 13485:2003
FM 40363

Ionto480 is a Registered Trademark of Axelgaard Manufacturing Co., Ltd.
05/10

©2010 Axelgaard Mfg. Co., Ltd.



Better Technology Better Results




Ionto480[®]

Buffered
Iontophoresis
Electrodes

Experience the Difference



AXELGAARD

MANUFACTURING CO., LTD.

Lonto480[®]

Buffered
Iontophoresis
Electrodes

Axelgaard Lonto480 Buffered Iontophoresis Electrode Treatment Kit is an innovative disposable iontophoretic drug delivery system, with Axelgaard's unique scientific blend of materials for optimum comfort during treatment.

- Each individually wrapped kit contains one active drug delivery electrode, one buffered return electrode and two alcohol swabs.

Active Drug Delivery Electrode

- Contains a reservoir pad to hold the deliverable solution.
- Reservoir pad fill area is conveniently accessible to allow for easy application of the solution prior to removal of the adhesive backing.
- The soft foam backing, with medical grade adhesive, is flexible and conforms to body contours.

Buffered Return Electrode

- The return electrode, with patented buffered gel, minimizes skin irritation.

Convenient Treatment System

*Lonto480
Active Drug
Delivery Electrodes
are Available in a Variety of Sizes*



Small Square

Drug Delivery Electrode
Approximate Fill Volume –
1.5cc

Medium Square

Drug Delivery Electrode
Approximate Fill Volume –
2.5cc

Large Square

Drug Delivery Electrode
Approximate Fill Volume –
4.0cc

Butterfly

Drug Delivery Electrode
Approximate Fill Volume –
2.0cc



**Return
Electrode**