


EconoStim®


- EconoStim® electrodes feature EconoStick® hydrogel, which is a thinner scrimless version of Axelgaard's world renowned patented MultiStick® hydrogel.
- The low-profile design of EconoStim enhances flexibility and allows the electrode to conform to all body contours.
- Innovative oversize adhesive border to maintain optimum adhesion and to facilitate ease of proper removal, which promotes electrode durability.
- In addition to the standard sizes, we have the EconoStim SELECT for TENS and stimulation of small muscle groups, microcurrent, and TMJ treatments.



EconoStim® Electrodes



Low Profile Electrode



Featuring EconoStick Hydrogel



Small Muscle Group NMES

These scrimless low-profile electrodes feature our patented EconoStick hydrogel. An excellent choice for TENS applications at an economical price.

Made in USA 



AXELGAARD
MANUFACTURING CO., LTD.


www.axelgaard.com
support@axelgaard.com
US + 1 (760) 451-8000
EU +45 8674 2088




Headquarters: Fallbrook, CA 92028 USA Europe: DK-8520 Lystrup, Denmark

Axelgaard Logo, MultiStick, EconoStim, and EconoStick are Registered Trademarks of Axelgaard Manufacturing Co., Ltd.
© 2019 Axelgaard Manufacturing Co., Ltd. 7/19

EconoStim®


available size and shape options



PART #	SHAPE	DESCRIPTION	#ELC/PK	#PK/CS	#ELC/CS
ESB5000		2.0" (5cm)	8	10	80
ESB5050		2.0" x 2.0" (5cm x 5cm)	8	10	80
ESB50100		2.0" x 4.0" (5cm x 10cm)	4	10	40



EconoStim®SELECT is well suited for TENS and stimulation of small muscle groups, microcurrent, and temporomandibular joint (TMJ) treatment.

PART #	SHAPE	DESCRIPTION	#ELC/PK	#PK/CS	#ELC/CS
ESB3050		1.2" x 2.0" (3cm x 5cm)	12	10	120

For further product information, please contact us.



AXELGAARD
MANUFACTURING CO., LTD.

www.axelgaard.com
support@axelgaard.com
US + 1 (760) 451-8000
EU +45 8674 2088

Headquarters: Fallbrook, CA 92028 USA Europe: DK-8520 Lystrup, Denmark

Axelgaard Logo, MultiStick, EconoStim, and EconoStick are Registered Trademarks of Axelgaard Manufacturing Co., Ltd.
© 2019 Axelgaard Manufacturing Co., Ltd. 7/19